

CORRIGENDUM 1

201.3.201

Replace the existing term " B_1 RMS" by " B_{1+} RMS".

Replace, in the definition:

*" B_1 , the radio frequency magnetic induction"
by*

*" B_{1+} , the MR relevant radiofrequency
magnetic induction"*

*Replace, in the formula, " B_1 RMS" by " B_{1+} RMS"
and " B_1 " by " B_{1+} ".*

201.3.201

Remplacer le terme existant " B_1 RMS" par " B_{1+} RMS".

Remplacer, dans la définition:

*" B_1 , l'induction magnétique à radiofréquence"
par*

*" B_{1+} , l'induction magnétique à radiofréquence
RM pertinente".*

Remplacer, dans la formule, " B_1 RMS" par " B_{1+} RMS" et " B_1 " par " B_{1+} ".

Table 201.101 – List of symbols

In the 2nd row, 1st column replace " B_1 " by " B_{1+} ".

In the 3rd row, 1st column replace " B_1 RMS" by " B_{1+} RMS".

In the 3rd row, 3rd column replace " B_1 " by " B_{1+} ".

Tableau 201.101 – Liste des symboles

Dans la 2^{ème} ligne, 1^{ère} colonne remplacer " B_1 " par " B_{1+} ".

Dans la 3^{ème} ligne, 1^{ère} colonne remplacer " B_1 RMS" par " B_{1+} RMS".

Dans la 3^{ème} ligne, 3^{ème} colonne remplacer " B_1 " par " B_{1+} ".

201.7.9.2.101 Instructions for use for MR EQUIPMENT

j) Exposure of the PATIENT to radio
frequency magnetic fields

In the final dashed item replace " B_1 RMS" by " B_{1+} RMS"

201.7.9.2.101 Instructions d'utilisation des APPAREILS A RM

j) Exposition du PATIENT aux champs
magnétiques à radiofréquence

Dans le dernier tiret, remplacer " B_1 RMS" par " B_{1+} RMS".